




# Domestic Scalable Plasma processing: technological approaches and solutions for LMIC

Josephine Cheng  
Senior Consultant Strategy Operationalization – traditional  
modalities – APAC, Merck Taiwan

Jian Fei Shou  
Senior MSAT Manager, Merck China

ISBT Workshop  
September, 2021



The life science business of Merck operates  
as MilliporeSigma in the U.S. and Canada.

## Millipore®

Preparation, Separation,  
Filtration & Testing Products 

## Disclaimer:

Views expressed in this talk constitute our professional opinion while being Merck Employees

**Millipore®**

Preparation, Separation,  
Filtration & Testing Products 

# Outline

1

**Overview on Plasma production**

2

**Safe plasma from safe processes- case study sharing**

3

**Proposals for strategic approach to improve domestic plasma production in LMIC**

4

**Take-home messages**

**Millipore®**

Preparation, Separation,  
Filtration & Testing Products



**Prosperous market outlook however unbalance supply for LMIC**

How to implement pragmatic and scalable technologies to make safe products to treat patients in need?

Wastage esp. in LMIC

Shortage

Intensified purification  
of Plasma IgG

How novel technologies  
helps to make safe  
PDMPs efficiently

Scalable process –every  
drop can be used

How to design a  
process from small to  
large scale

Process technology  
proposals

Single-use  
technologies accelerate  
flexible manufacturing

# From Plasma donation to patient administration Case study from public-private collaboration



MERCK

Blood Transfusion 2021;  
DOI 10.2450/2021.0159-21



## Merck White Paper

**Millipore**  
Purification, Separation  
Water & Wastewater Products

White Paper

### Intensification of Human Plasma IgG Purification for Intravenous and Subcutaneous Administration

Taipei Medical University, Taiwan:  
Thierry Burnouf, Distinguished Professor & Vice-Dean; Director, International Ph.D. Program in Biomedical Engineering; College of Biomedical Engineering  
Yu-Wen Wu, Post-doctoral Researcher, College of Biomedical Engineering

Merck:  
Josephine Cheng, Associate Director, Plasma & Vaccine Segment Development Asia  
Shiaren Wu, Manufacturing Science and Technology (MSAT) – Technology Specialist, Taiwan  
Karen Chan, Head of Downstream Manufacturing Science and Technology (MSAT), SEA and Taiwan  
Leo Liao, Senior Process Development Scientist, China  
Xisheng Cao, Senior Process Development Scientist, China  
Bin Wang, Head of Technical & Scientific Solutions, China

Therapeutic products derived from human plasma include coagulation factors, protease inhibitors, anticoagulants, albumin, polyvalent, and hyperimmune immunoglobulins (IgGs). IgGs are essential for treatment of patients with primary or secondary immune deficiencies and those with some inflammatory and autoimmune diseases.

The fractionation process used to extract proteins from plasma is designed to optimize recovery and ensure the appropriate quality and safety. In addition, the therapeutic IgG must be of a sufficient concentration for intravenous or subcutaneous administration and must meet stringent quality criteria including:

- Virus safety
- Low residual level of contamination by IgA, IgM, proteolytic enzymes, Factor XII/Xia or chemicals used for virus inactivation
- Lack of hemolytic effects due to the presence of anti-A and anti-B isoagglutinins

In this white paper, we present the results from our collaboration with Taipei Medical University to develop and evaluate the reliability and consistency of various new steps to purify plasma-derived IgGs. The intensified workflow includes flow-through mode chromatography and single-pass tangential flow filtration to achieve high recovery of a high quality product. The study focuses on the combination of purification steps to ensure a good removal of IgA, IgM, anti-A, anti-B, thrombogenic factors, and virus-inactivating agents.

The classical method of plasma purification is based on the cold ethanol precipitation approach developed decades ago by Edward Cohn. Today, fractionation is typically performed on a large scale with batch sizes of 500 – 10,000 liters of human plasma. Increasing concentrations of alcohol in the range of 8 – 40% are used to precipitate proteins according to their solubility at cold temperatures. The resulting Cohn fractions are crude starting materials requiring further purification based on a range of parameters including molecular size, charge, solubility and structure. IgG purification can start from either Fraction II, or with Fractions I, II and III to maximize the yield. In this case study, the intermediate IgG fraction used was purified from human plasma based on 20% pH 5.5 caprylic acid treatment, to represent a worse case scenario of Fraction I, II and III in terms of purity and relative proportion in IgG, IgA, and IgM.

The life science business of Merck operates as Merck KGaA in the U.S. and Canada.

Facilitating towards safe & efficient Plasma production

A generic, easy-to-operate, flowthrough-mode purification process that provides scalable & robust purification with enhanced productivity and quality IGG fitting for therapeutic usage.

### Quality criteria:

- Virus safety
- Low IgA & IgM contamination
- Low FXI/XIa
- Lack of Hemolytic effect
- Lack of chemicals used for virus inactivation


Millipore

## Summary key achievement of the collaboration




### Solid proof points for intensified IgG purification readily scalable for small-large-scale manufacturing





- ✓ FT-chrom for primary purification 82%→ 99%, mainly reduce IgA and IgM
- ✓ Eshmuno P reduced 8-32X anti-A/ anti-B isoagglutinins
- ✓ Triton-X 100/ TnBP provides strong Virus inactivation as soon as 5 minutes, and removal of S/D with Licroprep C18 resin in FT mode is effective (bdl)
- ✓ SPTFF technology (data not shown) provides a gentle way of inline/ final concentration, reaching 20% final target for SCIG purpose.
- ✓ All Aseptic filters/ prefilters showing robust filtration results and recovery (~ 100%)
- ✓ Overall process can be sliced & diced fitting the target end product(s), and can be easily incorporated a second virus removal method eg. Virus filtration, to meet regulatory requirement.
- ✓ All steps are readily scalable & implementable




# Proposal #1: Flowthrough one-step to remove IgA & IgM




Pilot scale in 2 batches  
IgG/IgA/IgM



Small scale in 10 cycles  
IgG1/IgG2/IgG3/IgG4



Average Subclasses (mean of 10 cycles)



## Summary #1 with AEX step:

1. Flow through one-step IgA/IgM removal
2. Purity IgG avg. 82% to 99% in small and pilot scale.
3. 200 cycles test
4. No changes in IgG subclasses.
5. No thrombogenicity activity detected


Learn More with our webinar:  
[Chromatography: Chromatographic strategies for IVIG purification - Part 2](#)




# Proposal #2:



## Robust reduction of the blood-type specific isoagglutinins




	Group A	Group B	Group AB	Group O
Red blood cell type				
Antibodies in Plasma			None	



**Eshmuno® P Anti-A (FT)**

**Eshmuno® P Anti-B (FT)**




- 8 to 16 times reduction in Anti-A titer
- 16 to 32 times reduction in Anti-B titer

\*samples tested at 30mg/ml concentration from 1<sup>st</sup> SPTFF (6X) step to the last step.



# Proposal #3: Second Virus inactivation step with S/D and removal by Licroprep C18 (40 – 63um)



## Key points:

- A. Classical TnBP/Triton X-100 provides > 4-5 LRV in time as short as 5 minutes.
- B. Typical chromatography for FT mode S/D-IGG running through C18 column, residual of S/D tested as low as 1ppm and 2ppm, respectively.

A.

Human IgG : 0.3% TnBP + 1% TX100 LRV Results					
Virus	Device	LRV at Incubation Time (min)			
		5	30	60	360
XMuLV	Mobius 1	≥5.5	≥5.3	≥5.3	≥5.4
	Mobius 2	≥5.5	≥5.3	≥5.3	≥5.5
BVDV	Mobius 1	≥4.5	≥4.4	≥4.6	≥4.5
	Mobius 2	≥4.4	≥4.6	≥4.4	≥4.5

B.

Residual TnBP of SD-IgG (Ratio of resin and loaded IgG)	Batch 3 (1mL:6mL)
<b>C18</b>	<1 ppm
<b>SPTFF-5X</b>	<1 ppm

Residual Triton X-100 of SD-IgG (Ratio of resin and loaded IgG)	Batch 3 (1mL:6mL)
<b>C18</b>	<2 ppm
<b>SPTFF-5X</b>	<2 ppm

**Learn More with our webinar:**  
[Solvent Detergent Viral Inactivation using S.U Technology in Blood Fractionation Processes](#)

\*Table A source: Hsieh YT, Mullin L, Greenhalgh P, Cunningham M, Goodrich E, et al.: Single-use technology for solvent/detergent virus inactivation of industrial plasma products. *Transfusion* 2016; 56: 1384-93.

# Outline

1

**Overview on Plasma production**

2

**Safe plasma from safe processes- case study sharing**

3

**Proposals for strategic approach to improve domestic plasma production in LMIC**

4

**Take-home messages**

**Millipore®**

Preparation, Separation,  
Filtration & Testing Products



A pragmatic approach

Establishing domestic capability on plasma processing

MERCK

Understand cost structure

Think scalable from the start


Incorporate single-use to facilitate competency

Millipore 



# Cost Analysis



# Understanding the cost structure in IgG processing





200L Hybrid Cost Breakdown



Capital	11.28 USD Million
Cost of Goods	108.1 USD/g
PMI	8,165
Capacity	198.0 kg/yr
Doses per Year	7.9E+04


2000L Hybrid Cost Breakdown





Capital	24.21 USD Million
Cost of Goods	18.3 USD/g
PMI	3,905
Capacity	1979.9 kg/yr
Doses per Year	7.9E+05

- Costly steps are IEX, AC, VF, and UF.
- In a hybrid (SU-SS) process, labor cost ~20-30% depending on scale.

# Understanding the cost structure in IgG processing





200L Stainless Steel Cost Breakdown



Capital	25.91 USD Million
Cost of Goods	157.6 USD/g
PMI	33,796
Capacity	191.0 kg/yr
Doses per Year	7.6E+04

2000L Stainless Steel Cost Breakdown




Capital	34.77 USD Million
Cost of Goods	21.9 USD/g
PMI	8,991
Capacity	1910.4 kg/yr
Doses per Year	7.6E+05

- Costly steps are IEX, AC, VF, and UF.
- In a non-SU process, labor cost ~40-50% depending on scale.



### Cost per dose at different production scale (USD/g) reduces when scale increases



- Unit operations in a IGG purification process which require the most cost are IEX, AC, VF, and UF. Optimizing these steps can make the most impact on reducing the overall cost.
- Incorporating Single-use technologies reduce the Capex and Labour cost, mainly due to the elimination of large systems, and cleaning/validation time.
- Though consumable cost will be higher, the overall cost of including SU technologies can help to accelerate time of establishment, time to train employee, reduced footprint needed, and eliminate risk of human error related contamination.



## WHY Single-Use?



## Examples in single-use

### Mixing



### Assemblies, Connectors, Samplings



### Storage & Transport



### Final Filling



# Imaging a site with simplicity



Stainless Steel

VS

Single-Use



# Flexible & Next generation manufacturing taking vaccine as example

## Implementation of Single use in Final filling – GSK case study H1N1, 2009


**Faster deployment**  
**Flexibility to change scale or process**  
**Reduces time to market**  
**Accelerates response to high surge of vaccines: this can well apply for PLASMA**

	Traditional	Single use
Clean and set-up	14 Hrs	<1 Hr
Cleaning validation	Extensive	Zero
Filling time	24 hrs	10 hrs
Average vials/hr	3,000	10,000
Aseptic connections	50	0
Operator training	2 weeks	2 days
Equipment utilization	35%	82%
<b>Total time</b>	<b>38 hrs</b>	<b>12 hrs</b>

Traditional large vaccine manufacturing facilities



Manufacturing facility using single-use technologies



	Traditional stainless facility	Single-use facility
Capex required	~\$500M to \$1B	\$20-100M
Time to construct	5-10 years	1.5 years
Change over time	4 weeks	0.5 days
Footprint	~>70,000 m <sup>2</sup>	~11,000 m <sup>2</sup>

Scalable process eliminates all kinds of  
process upsets and surprises





# Scalability









Figure 5. Throughput scalability for Vis



- Mini Avg
- Mini Avg+10%
- Mini Avg - 10%
- C5EN89074-3137
- C5EN89074-3139
- C5EN89074-3152
- C5EN89076-2912
- C5EN89076-2932
- C5EN89076-2943
- C5EN89080-2644
- C5EN89079-2722
- C5EN89079-2726
- C5EN89079-2730
- C5EN89080-2608
- C5EN89080-2616
- C6AN44268-5075
- C6AN44268-5066
- C6AN44268-5073
- C6BN44270-5109
- C6BN44270-5116
- C6BN44270-5123
- C6BN44271-5146
- C6BN44271-5150
- C6BN44271-5159
- C6BN44272-5170
- C6BN44272-5172
- C6BN44272-5177
- 33 m²
- 1.1 m²
- 23 cm², #1
- 23 cm², #2
- 23 cm², #3
- 23 cm², #4
- 540 cm², #1
- 540 cm², #2
- 540 cm², #3
- 540 cm², #4






Does Single-use scale up well?




# Performance comparison

## Pellicon capsule (single-use) vs cassette (multi-use)


Comparable flux and pressure drop




● Capsule 0.1 m²    ● Capsule 0.5 m²  
 □ Cassette 0.11 m²



Linear scalability within capsule family and to cassettes




[SPTFF] Comparable single-pass conversion



# How We Serve the Plasma Industry


**MERCK**

Products  
Quality dossiers




Process Intensification  
& downstream efficiency

Services  
Process  
development



Speed to clinic and simplified compliance

Collaborations  
Training



Industry and academic  
partnerships

**Millipore** 



1

Increasing self-sufficiency for PDMPs to address the global demand is a must, domestic production is a key starting point.

2

Product development can start simple (e.g. S/D treated cryo) to more complex (e.g. IgG, Albumin, Factor VIII).  
Think scalable when developing a process to ensure smooth progression from small (e.g. 100L) to large scale (e.g. 2000L); incorporate disposable technologies to accelerate fractionation competency establishment.

3

Collaboration accelerates development of plasma production even starting from small scale



# Acknowledgements

## Taipei Medical University, Taiwan (TMU):

- Prof. Thierry Burnouf
- Yu-Wen Wu
- Chen-Yun Wang
- Cheum Lam Hong

## Merck:

- Sharon ShangJung Wu
- Karen Waiyu Chan
- Leo Xun Liao
- Xisheng Cao
- Bin Wang
- Ashok Kumar

**Millipore®**

Preparation, Separation,  
Filtration & Testing Products



# Thank You

# Q&A

[Josephine.cheng@merckgroup.com](mailto:Josephine.cheng@merckgroup.com)

[Jian.fei.shou@merckgroup.com](mailto:Jian.fei.shou@merckgroup.com)



Connect to speaker



[http://bit.ly/26SM\\_LNJC](http://bit.ly/26SM_LNJC)



Connect to speaker

