



Potential Contribution of Nanofiltration to the Virus Safety of Domestic Plasma Components

September 23 2021

Tomoko Hongo-Hirasaki,

Lead Expert, Bioprocess Division, Asahi Kasei Medical Co., Ltd.

CONFIDENTIAL © 2021 Asahi Kasei Bioprocess



Tomoko Hongo-Hirasaki is an employee: a lead expert of Bioprocess Division, Asahi Kasei Medical Co., Ltd., Tokyo, Japan.



Purpose of this presentation

- To provide an understanding of nanofiltration
- To reduce barriers to the use of nanofiltration
- To facilitate the use of nanofiltration as a virus removal process



1. Nanofiltration
2. Challenges for Implementation
3. Case Study of Implementation
4. Conclusion



1

Nanofiltration



- The most robust virus removal methods recognized by the plasma fractionation industry and in international regulatory guidance.
- Applicable to various types of plasma products.
- Can remove various known and emerging viruses based on a size-exclusion mechanism.
- Excellent scalability: straightforward process for implementation from mini-pools to large-scale batches of plasma fractions.



- The most robust virus removal methods are recognized by the plasma fractionation industry and in international regulatory guidance.
- Applicable to various types of plasma products
- Can remove various known and emerging viruses based on a size exclusion mechanism.
- Excellent scalability: straightforward process for implementation from mini-pools to large-scale batches of plasma fractions.



- Viral safety of bio-therapeutic products is required by regulatory agencies worldwide.
- Viral inactivation/ removal method(s) are necessary in production processes.
- “It is desirable to investigate the contribution of more than one production step for virus reduction and at least two orthogonal steps should be assessed”*

* EMEA/ CHMP/ BWP/ 398498, 2005



- Nanofiltration is a robust virus removal method

WHO Technical Report, Series No. 924, 2004

Annex 4

Guidelines on viral inactivation and removal procedures intended to assure the viral safety of human blood plasma products



The European Agency for the Evaluation of Medicinal Products
Evaluation of Medicines for Human Use

London, 28 March 2001
CPMP/BWP/BPWG/4080/00

EMEA WORKSHOP ON VIRAL SAFETY OF PLASMA-DERIVED MEDICINAL PRODUCTS WITH PARTICULAR FOCUS ON NON-ENVELOPED VIRUSES

13 September 2000
REPORT

Received: 20 November 2019 | Revised: 24 June 2020 | Accepted: 25 June 2020

DOI: 10.1111/trf.16022

BLOOD COMPONENTS

TRANSFUSION

Nanofiltration as a robust method contributing to viral safety of plasma-derived therapeutics: 20 years' experience of the plasma protein manufacturers

Nathan J. Roth¹ | Herbert O. Dichtelmüller^{2†} | Fabrizio Fabbrizzi^{3†} |
Eckhard Flechsig² | Albrecht Gröner⁴ | Mary Gustafson⁵ | Juan I. Jorquera⁶ |
Thomas R. Kreil⁷ | Dominika Misztela⁸ | Elisa Moretti³ | Mila Moscardini³ |
Gerhard Poelsler² | John More⁹ | Peter Roberts⁹ | Andreas Wieser⁷ |
Rodrigo Gajardo⁶



- The most robust virus removal methods are recognized by the plasma fractionation industry and in international regulatory guidance.
- Applicable to various types of plasma products.
- Can remove various known and emerging viruses based on a size-exclusion mechanism.
- Excellent scalability: straightforward process for implementation from mini-pools to large-scale batches of plasma fractions.





Planova filter



Hollow fiber




Membrane cross-section



- Hollow fibers of Planova filters have a three-dimensional network structure consisting of voids connected by capillaries.




Filtration Mechanism: Size Exclusion



Planova Filters Against Emerging Viruses

- Virus removal based on size exclusion mechanism



- Planova filters are applicable to a wide range of plasma products



Pharmaceuticals	M.W. (kD)	15N	20N/ BioEX	35N
Thrombin	30 ~ 40			
α 1Protease inh.	56			
Factor-IX	58			
Protein C	62			
Hemoglobin, antithrombin	65			
Albumin	66			
polyclonal IgG	160			
Factor XI	180			
Factor-VIII	300			
Fibrinogen	340			
IgM	800			



- The most robust virus removal methods are recognized by the plasma fractionation industry and in international regulatory guidance.
- Applicable to various types of plasma products.
- Can remove various known and emerging viruses based on a size exclusion mechanism.
- **Excellent scalability: straightforward process for implementation from mini-pools to large-scale batches of plasma fractions.**



- Viral clearance and flow rate scalability data for Planova filters



Brian Buesing, Asahi Kasei Bioprocess America, *Planova Workshop*, Athens, 2015 (excerpted)

- Robust viral clearance across both filter types and sizes
- Excellent filter: lot-to-lot consistency in virus removal
- The linear increase in flow rates shows uniform scalability of flux




2

Challenges for Implementation



Main Challenges

- ✓ Introducing nanofiltration in a process:
 - How to use it?
 - How to optimize the filtration conditions?
- ✓ How to meet regulatory requirements and technological trends regarding virus safety






Proposed Solutions

- ✓ Hands-on training to use Planova filter
- ✓ Providing in-depth technical advice and support services
- ✓ Providing regulatory information and updated scientific knowledge sharing
- ✓ On-site training on virus clearance study




Lab test



Non-GLP

GLP (Validation)

Virus clearance study



Regulatory approval

Pilot test

Commercial production

Scale-up



- Annual event where “Planova Family” gathers to exchange information
- Held since 1998, alternately in Europe and the United States
- Planova users, regulatory authorities and key opinion leaders present their experience of the regulation, downstream process development, pathogen safety and the manufacturing of biopharmaceuticals
- Asahi provides new data or other technical updates on nanofiltration



Planova Workshop, Rome, 2013



- Planova filters have been used for biological products for almost 30 years worldwide, including LMIC



3

Case Study of Implementation



Optimal operating filtration conditions

Filter size: 0.001m²
Operating Pressure, temperature, protein concentration,
Prefilter

Process scale up and product integrity

Filter size: 0.1m²
Operating Pressure: 294 kPa
Temperature: 25°C
Protein concentration: 50 g/L
Prefilter: 0.1 μm


Preliminary viral clearance study

Filter size: 0.001m²

Viral clearance validation study at NIFDC

Filter size: 0.001m²


- This plasma derived IVIG product is the first to be certified for parvovirus safety by the NIFDC in China.




Sha Ma et.al., Biologicals 52(2018) 37-43 (adapted)




Flow chart of preliminary test: Experimental design



CA-IgG nanofiltration by Planova 35N



CA-IgG nanofiltration by Planova 20N



- ✓ Planova 35N nanofiltration:
50-60 mL of CA-IgG could be readily nanofiltered on 0.001 m² Planova 35N within 2hrs at a constant pressure of 0.8 kgf/cm²
- ✓ Planova 20N nanofiltration:
16 mL of CA-IgG could be nanofiltered on 0.001m² Planova 20N within 2hrs at a constant pressure of 0.8 kgf/cm²
- ✓ Since CA is a virus inactivation step (enveloped viruses), nanofiltration of CA-IgG provides a second virus reduction step (enveloped and non-enveloped viruses) to consider when batch size increases



4

Conclusion



- Implementing nanofiltration is a straightforward process from mini-/small-pools domestic plasma components (e.g., immunoglobulins) to large-scale batches of plasma fractions (e.g., Factor VIII, Factor IX, immunoglobulins).
- Technical support and experienced/scientific knowledge sharing can be provided to LMIC to improve the virus safety of domestic plasma components.



AsahiKASEI

BIOPROCESS

감사합니다

Gracias

Danke

Благодаря

谢谢

Tack

धन्यवाद

Dziękuję

Спасибо

Thank You

Obrigado

Děkuju

Grazie

Ευχαριστώ

Merci

Köszönöm

ありがとうございました

Teşekkür ederim

